

Print & Share Journal

The Report screen lets you generate a summary of your readings, meal log, or medication notes. Share with anyone you like — print a paper copy, save as PDF, or email directly.

Step-by-Step: Generate a Report

1

Open Report Screen

From the Home screen, tap the Print / Share journal.

2

Choose Log Type

Select Sugar Log (readings), Food Log (meals), or Medication. You can select multiple log types to share together.

3

Set Date Range

Tap the From and To fields to set your range.

4

Share

You can print, share as PDF, or email your contacts.

5

Preview

Tap Preview Journal to see what your report will look like before sharing.

■ If My Contacts is filled out with email addresses, these will be saved to your device clipboard for easy copy/paste inside your email app.

GD SugarLog is a personal tracking tool and is not a substitute for medical advice.